

Process-suited Raman Spectroscopy with High Optic Collection Efficiency for In-line Analytics

F. Braun¹, S. Schwolow², R. Schalk¹, J. Seltenreich¹, M. Nachtmann¹, D. Geörg¹, T. Röder², T. Beuermann¹, N. Gretz³, M. Rädle¹

¹Hochschule Mannheim, Institute for Process Measurement Technology and Innovative Energy Systems, Mannheim

²Hochschule Mannheim, Institute for Chemical Process Engineering, Mannheim

³Heidelberg University, Center for Medical Research, Medical Faculty Mannheim

INTRODUCTION

Established off-line methods for industrial monitoring of chemical reactions increasingly replaced by optical in-line methods


- ➔ Application of in-line Raman spectroscopy with fiber-optic probe
- ➔ Significant time and material savings compared to off-line analytics (eliminating additional work steps)
- ➔ No safety or contamination risks during handling

Spectrometer optics optimized for maximum photon flux

- ➔ Rapid sampling rates despite reduced laser output


Michael addition as a quick sample reaction with known kinetics [1]

- ➔ Rapid recording of time conversion curves in the batch reactor
- ➔ Low excitation power (7 mW)




MEASUREMENT TECHNOLOGY

Backscatter probe, immersed (no focusing optics)



Raman spectrometer optics (custom-built) suited to the application


- ➔ **Laser light source**
Excitation wavelength: 532 nm -> limitation: fluorescence background
- ➔ **Spectrometer resolution**
Variable through the use of different detection fibers
- ➔ **Process suitability**
Robust; spectrometer without movable components and no filters, mirrors, or lenses in the probe



RESULTS


Measurements in the batch reactor under backscattering

- ➔ Excitation power: 7 mW; integration time: 1000 ms; accumulation: 1
- ➔ Evaluation by vinyl band at 1637 cm⁻¹



Monitoring a Michael addition reaction (n=5)


- ➔ Raman measurements and comparison with kinetic model



- ➔ Good agreement with kinetic model based on tests with GC off-line analytics [2]

CONCLUSION

- ➔ Realization of spectrometer optics optimized for maximum photon flux
- ➔ Efficient methods for rapid measurement of kinetics through in-line reaction monitoring
- ➔ Kinetics of model reaction monitoring confirmed by conventional Raman measurement technique [1] and GC analysis [2]
- ➔ Allowable continuous power output for laser light sources in explosion-protected areas as per ATEX (DIN EN 60079-28) of 7 mW is below limit by a factor of 5 -> great potential in operational practice
- ➔ Range of applications can markedly expand in conjunction with chemometric methods



Literature

- [1] Schwolow, S., Braun, F., Rädle, M., Kockmann, N., Röder, T., 2015. Fast and Efficient Acquisition of Kinetic Data in Microreactors Using In-Line Raman Analysis. Org. Process Res. Dev.
- [2] Schwolow, S., Heikenwälder, B., Abahmane, L., Kockmann, N., Röder, T., 2014. Kinetic and Scale-up Investigations of a Michael Addition in Microreactors. Org. Process Res. Dev.